

## APPLICATION NOTE

### Validation of a robust PCR-only assay for quantifying *FMR1* CGG repeats on clinical samples

#### Authors

**J. Kent Moore**  
Galina Filippov  
Karl Adler  
Mack Schermer  
PerkinElmer, Waltham, MA, USA

**Monica J. Basehore**  
Michael J. Friez  
Greenwood Genetic Center  
Greenwood, South Carolina, USA

**Heike Fiegler**  
PerkinElmer, UK

## FragilEase - A New PCR-based Assay For Fragile X Testing

### Introduction

Fragile X syndrome (FXS) is the most common form of inherited intellectual disability and is caused by an expansion of a CGG repeat in the 5' end of the *FMR1* gene. Fragile X testing is traditionally performed using a *FMR1*-specific PCR

followed by electrophoresis using a capillary sequencer. Typical laboratory developed PCR approaches have been hampered however by the ability to reliably amplify the CGG-rich region beyond about 100 – 130 repeats. Furthermore, electrophoresis with a standard capillary sequencer is capable of measuring *FMR1* PCR products only up to approximately 200 repeats with commercially available size standards.


An accurate measurement of the repeat size however is required for correct classification and diagnosis of FXS and other Fragile X-associated disorders. Thus, differentiating full mutations with greater than 200 CGG repeats from homozygous normal female samples, and confirming full mutations, has historically required labor-intensive measurement by Southern blot.

For this study we have used the FragilEase™ kit, a simple and robust PCR-based assay for quantification of CGG repeat alleles together with a low-cost benchtop electrophoretic instrument. The FragilEase™ PCR assay is designed to amplify the entire CGG repeat sequence in the *FMR1* promoter region thus allowing an accurate amplification and reliable detection of normal, intermediate, pre- and full mutations.

Not for distribution in the USA.

For Research Use only. Not for Diagnostics use.  
Please check product availability with your PerkinElmer sales representative.

## Workflow



### Fragile X testing using FragilEase™

A blinded set of 200 archived DNA samples was analyzed by use of the FragilEase™ kit as previously described (1). The samples had previously been characterized by standard laboratory PCR and reflex testing by Southern analysis (2,3). The sample set included both male and female individuals with normal (n=167), intermediate (n=10), premutation (n=11), and full mutation alleles (n=12).

After PCR completion, the PCR fragments were purified using the PureLink Micro PCR purification kit (Invitrogen, USA). The purified PCR fragments were then analyzed using a 2100 Bioanalyzer (Agilent, USA) instrument. Two Coriell reference samples (Coriell, USA) with well characterized repeat sizes were included as standards in each set-up.

### Analysis using FraXsoft™

Using FraXsoft™, the repeat number for each sample was calculated based on the measured fragment size. Including reference samples with well characterized repeat sizes in each experiment allows the construction of a linear regression standard curve. The measured fragment sizes of all test samples are then plotted against this linear regression curve, allowing an accurate determination of the repeat size of each test sample. According to the characterization guidelines chosen, the mapped repeat size will be displayed in its corresponding categories (i.e. normal, intermediate, premutation, full mutation).




Figure 1. Two example cases representing representing full mutations samples are shown for a typical Fragile X Southern Blot together with the Bioanalyzer electropherogram and FraXsoft report after DNA amplification with FragilEase™.

## Results

Of the 200 DNA samples tested, 199 samples were concordant with the previously determined classification. All twelve full-mutation samples were successfully classified, with the largest allele size measured at over 653 repeats. The first round PCR failure rate was 2% (4/200) with all four samples producing correct classification result on a subsequent second round of PCR.

The only one discordant sample was previously identified by Southern blot analysis to have 45 repeats and was subsequently classified as intermediate according to the

ACMG guidelines. In this study, the calculated repeat number for the same sample was found to be 44 repeats. Using the ACMG guidelines this would classify as normal (normal <45 repeats). However, the precision for allele size detection using FragilEase™ together with the Bioanalyzer ranged from +/-1 to +/-2 repeats for alleles with repeat sizes of up to 58 repeats in our study. Hence, this sample is within the detection limits and we would recommend confirmation by an independent methodology in similar cases.

### Male Sample Results

New FragilEase™ PCR-only Method						Reference Method
Classification	Normal	Intermediate	Premutation	Full Mutation	Total	
Normal	87	0	0	0	87	
Intermediate	1	6	0	0	7	
Premutation	0	0	0	0	0	
Full Mutation	0	0	0	6	6	
<b>Total</b>	<b>88</b>	<b>6</b>	<b>0</b>	<b>6</b>	<b>100</b>	

### Female Sample Results

New FragilEase™ PCR-only Method						Reference Method
Classification	Normal	Intermediate	Premutation	Full Mutation	Total	
Normal	80	0	0	0	80	
Intermediate	0	3	0	0	3	
Premutation	0	0	11	0	11	
Full Mutation	0	0	0	6	6	
<b>Total</b>	<b>80</b>	<b>3</b>	<b>11</b>	<b>6</b>	<b>100</b>	

## Summary

The FragilEase™ assay was used on 200 blinded samples enriched for a majority of normal samples. All samples were correctly classified including eleven premutation and twelve full mutation samples. Accurate measurement of the CGG trinucleotide repeat in the *FMR1* gene is part of the classification and diagnosis of the Fragile X syndrome and other Fragile X-associated disorders. While this assay doesn't determine methylation status of the patient sample, the ability to robustly amplify premutations and more importantly, full mutations, will reduce the requirement of reflex testing by Southern analysis. Furthermore, this assay fulfills an

unmet need for a quick, low-labor and low equipment cost assay, allowing Fragile X testing to be performed in various laboratories within one working day.

## References:

- [1] Adler et al., 2011
- [2] Basehore et al., 2012
- [3] Basehore and Friez, 2009

Not for distribution in the USA.

For Research Use only. Not for Diagnostics use.

Please check product availability with your PerkinElmer sales representative.

PerkinElmer, Inc.  
940 Winter Street  
Waltham, MA 02451 USA  
Phone: (800) 762-4000 or  
(+1) 203-925-4602  
www.perkinelmer.com

PerkinElmer, Inc.  
Wallac Oy  
PO Box 10  
20101 Turku, Finland  
Phone: (+358) 22678-111  
Fax: (+358) 22678-357

ISO 13485  
ISO 9001  
CMDCAS  
ISO 14001  
OHSAS 18001



For a complete listing of our global offices, visit [www.perkinelmer.com/ContactUs](http://www.perkinelmer.com/ContactUs)

©2013 PerkinElmer, Inc. All rights reserved. PerkinElmer is a registered trademark of PerkinElmer, Inc. All other trademarks depicted are the property of their respective holders or owners. PerkinElmer reserves the right to change this document at any time and disclaims liability for editorial, pictorial or typographical errors.

1599-9776, December 2013